



# Sensitive, Quick LC/MS/MS Analysis of Aflatoxins in Cannabis


## Using Supel<sup>™</sup> Tox AflaZea SPE and an Ascentis<sup>®</sup> Express Phenyl-Hexyl HPLC Column


Quantitative analysis of mycotoxins commonly involves sampling, sample preparation, extraction, and cleanup followed by chromatographic methods such as GC and HPLC. The matrix complexity of cannabis often makes sample cleanup methods used for common commodities ineffective. In this study, a sample preparation method using Supel™ Tox AflaZea SPE cartridges for cleanup of cannabis extracts prior to LC/MS/MS analysis of aflatoxins was developed. Supel™ Tox AflaZea cartridges utilize the "interference removal" strategy, requiring few processing steps and saving time by eliminating wash steps prior to analyte elution. The Ascentis® Express Phenyl-Hexyl HPLC column provided the selectivity to separate four aflatoxin compounds in cannabis.

### **Experimental**

Dried cannabis sample was obtained courtesy of Dr. Hari H. Singh, Program Director at the Chemistry and Physiological Systems Research Branch of the National Institute on Drug Abuse at the National Institute of Health. The sample was ground to a fine powder and extracted following the procedure outlined in the condition section of **Figure 2**. Pictures of the cannabis samples before, during, and after cleanup are shown in **Figure 1**. Matrix-matched calibration curves were constructed and run along with solvent-based calibration curves to compare ionization effects and sample cleanliness.

Figure 1. Photos of the Cannabis Samples (A) Before Cleanup, (B) On Supel™ Tox AflaZea SPE Cartridge, and (C) After Cleanup









#### **Results and Discussion**

#### **Matrix Removal**

The chromatographic separation of the aflatoxins was performed on an Ascentis® Express Phenyl-Hexyl HPLC column (**Figure 2**). Cannabis matrix effects were determined by comparison of the calibration curves constructed in solvent versus those in extract. Significant ion suppression was observed in the cannabis samples, and the matrix-matched calibration curves were required for accurate quantitation. The matrix effects can be attributed to the complex cannabis composition and the limited capacity of the Supel™ Tox AflaZea SPE for removal of all of the components.

Analyte recovery values from spiked cannabis extracts fall in the range of 102-127% with RSD below 12% for three replicates (**Table 1**). Excellent recovery values were even observed for Aflatoxins  $B_2$  and  $G_2$  that were spiked at a low level of 6.1 ppb.

Table 1. Percent Recovery for Aflatoxins from Cannabis (n=3)

	Aflatoxin B <sub>1</sub>	Aflatoxin B <sub>2</sub>	Aflatoxin G <sub>1</sub>	Aflatoxin G₂
Recovery (%)	102	109	108	127
RSD% (n=3)	8	12	3	9

<sup>\*</sup> versus matrix-matched calibration curve

#### **Conclusions**

A sample preparation method utilizing Supel™ Tox AflaZea SPE cartridges was developed for the cleanup of cannabis extracts prior to HPLC analysis. SPE, used according to the standard methodology, contributed to a simple, economical, quick analysis. This SPE methodology, in combination with LC/MS/MS detection and the selectivity of the Ascentis® Express Phenyl-Hexyl HPLC column, allowed for sensitive detection of four aflatoxin compounds in cannabis with recoveries of 102−127% at 6−25 ppb levels. Therefore, this current analytical method utilizing SPE, UHPLC, and MS/MS detection can be used successfully for testing aflatoxins at 5–50 ppb in cannabis.

#### **Featured and Related Products**

Description	Cat. No.
SPE Tube	
Supel™ Tox AflaZea SPE Cartridge 6 mL, pack of 30	55314-U
HPLC Column	
Ascentis® Express Phenyl-Hexyl 5 cm $\times$ 2.1 mm I.D., 2.7 $\mu$ m particle size	53334-U
Standard	
Aflatoxin Mix 4 Solution 0.5 $\mu$ g/mL B <sub>2</sub> and G <sub>2</sub> in acetonitrile 2 $\mu$ g/mL B <sub>1</sub> and G <sub>1</sub> in acetonitrile	34036
Solvent	
Acetonitrile, for HPLC, ≥99.9%	34851

Visit our cannabis testing resources at **SigmaAldrich.com/cannabis** 

#### Distributed By:

Greyhound Chromatography and Allied Chemicals 6 Kelvin Park, Birkenhead, Merseyside, CH41 1LT, UK

Tel: +44 (0)151 649 4000 Fax: +44 (0)151 649 4001

Email: sales@greyhoundchrom.com

WWW.GREYHOUNDCHROM.COM





