


MACHEREY-NAGEL

NUCLEODUR[®] π²

Chromatography



Now available in 5 μm and 3 μm particle size

- Hydrophobic phase with alternative selectivity compared to classical C₁₈ modifications
- Enhanced π-π interactions in HPLC
- Suitable for 100 % aqueous mobile phase



MACHEREY-NAGEL

www.mn-net.com



Key features

- Hydrophobic phase with alternative selectivity compared to classical C₁₈ modifications
- Separation mode based on two retention mechanisms: π - π interactions and hydrophobic interactions
- Better retention of aromatic and unsaturated substances
- Excellent performance under highly aqueous conditions
- Pronounced orthogonal selectivity in acetonitrile and especially in methanol

Recommended applications

- Overall sophisticated analytical separations
- Aromatic and unsaturated compounds
- Polar compounds
- Pharmaceuticals
- Antibiotics
- Steroids
- Pesticides



USP L11

Similar phases: Pinnacle® DB Biphenyl, Ultra Biphenyl

Technical data

Biphenylpropyl modification with multi-encapping

pH stability:	3–10
Particle size:	3 μm and 5 μm
Pore size:	110 Å
Specific surface:	340 m ² /g
Carbon content:	17 %



Stability under aqueous conditions

MN Appl. No. 127900

Chromatographic conditions


Column:	125/4 NUCLEODUR® π^2 , 5 μm
MN REF:	760623.40
Eluent:	20 mmol/L potassium dihydrogen phosphate in water (adjusted to pH 3.1 by 85 % phosphoric acid)
Flow rate:	1 mL/min
Temperature:	20 °C
Detection:	UV, 220 nm
Injection:	1 μL cytosine in methanol–buffer (1:1, v/v)

Enhanced stability

Excellent performance using 100 % aqueous mobile phase.




Chromatogram



Steroids

MN Appl. No. 127910

Chromatographic conditions

 Column: 125/4 each
 NUCLEODUR® π², 5 μm
 NUCLEODUR® Phenyl-Hexyl, 5 μm
 NUCLEODUR® C₁₈ Gravity, 5 μm
 MN REF: 760623.40 (NUCLEODUR® π²)
 Eluent: acetonitrile–water (45:55, v/v)
 Flow rate: 1 mL/min
 Temperature: 25 °C
 Detection: UV, 230 nm
 Injection: 1 μL


Analyte	Peak No.
1. Estriol	1
2. Hydrocortisone	2
3. Prednisone	3
4. β-Estradiol	4
5. Corticosterone	5
6. Cortisonacetate	6
7. Testosterone	7

Alternative selectivity

A complete separation was only achieved on NUCLEODUR® π².




Chromatograms



Amines

MN Appl. No. 128580

Chromatographic conditions

 Column: EC 150/2 NUCLEODUR® π², 5 μm
 MN REF: 760624.20
 Eluent A: 0.02 % formic acid in water
 Eluent B: 0.02 % formic acid in acetonitrile
 Gradient: in 10 min from 5 % to 85 % B, hold for 10.0 min,
 in 1.0 min to 5 % B, hold 5 % B for 9.0 min
 Flow rate: 0.3 mL/min
 Temperature: 40 °C
 Detection: MS, SMRM
 Injection: 5 μL


Concentration

300 ng/mL for each analyte

MRM transitions

Analyte	RT [min]	[M+H] ⁺	Q ₁ (Quantifier)	Q ₂ (Qualifier)
Trimethylamine	0.92	60.0	44.0	45.0
2-(Diethyl-amino)-ethanol	0.96	118.0	100.1	72.0
Pyridine	1.08	80.0	53.0	
Aniline	1.49	93.9	77.1	
p-Toluidine	2.11	107.9	93.1	65.1

Chromatogram




Analyte	RT [min]	[M+H] ⁺	Q ₁ (Quantifier)	Q ₂ (Qualifier)
N,N-Dimethyl-benzylamine	2.22	136.0	91.0	65.0
Diisobutylamine	2.37	136.0	91.0	65.0
Dibutylamine	3.09	130.1	57.0	41.0
Dibenzylamine	8.87	198.0	91.1	65.0

Anthocyanins in blueberry juice

MN Appl. No. 128260


Chromatographic conditions

	Column:	EC 250/4 NUCLEODUR® π ² , 5 μm
	MN REF:	760625.40
	Eluent A:	5 % formic acid in water
	Eluent B:	methanol
	Gradient:	5–95 % B in 30 min, hold for 10 min, back to 5 % B in 1 min, hold for 9 min
	Flow rate:	0.75 mL/min
	Temperature:	65 °C
	Detection:	Vis, 520 nm
	Injection:	20 μL

Concentration

Blueberry juice was diluted 10 times with eluent A, following filtration with CHROMAFIL® Xtra PET-20/13 (MN REF 729222)

Chromatogram



MRM transitions


Analyte	RT [min]
Delphinidin 3-galactoside	10.5
Delphinidin 3-glucoside	10.8
Delphinidin 3-arabinoside	11.6
Cyanidin 3-galactoside	11.9
Cyanidin 3-glucoside	12.2
Petunidin 3-galactoside	12.5
Petunidin 3-glucoside	12.7

Analyte	RT [min]
Cyanidin 3-arabinoside	13.0
Petunidin 3-arabinoside	13.6
Peonidin 3-galactoside	13.9
Malvidin 3-galactoside + Peonidin 3-glucoside	14.2
Malvidin 3-glucoside	14.4
Malvidin 3-arabinoside	14.9

Triarylmethane dyes

MN Appl. No. 128430


Chromatographic conditions

	Column:	EC 100/3 NUCLEODUR® π ² , 3 μm
	MN REF:	760636.30
	Eluent A:	5 mmol/L ammonium acetate + 0.1 % formic acid in water
	Eluent B:	acetonitrile
	Gradient:	in 10 min from 20 % to 85 % B, hold for 5.0 min, in 1.0 min to 20 % B, hold 20 % B for 5.0 min
	Flow rate:	0.4 mL/min
	Temperature:	25 °C
	Detection:	MS, MRM
	Injection:	5 μL

Concentration

0.25 ng/mL in eluent A for each analyte

Chromatogram



MRM transitions


Analyte	RT [min]	[M+H] ⁺	Q ₁ (Quantifier)	Q ₂ (Qualifier)
Malachite Green	10.1	328.9	313.1	207.9
Malachite Green-d ₅	10.1	333.9	318.3	213.1
Leucocystal Violet	11.2	373.9	358.0	238.2
Leucomalachite Green-d ₅	11.7	335.9	321.1	243.1

Analyte	RT [min]	[M+H] ⁺	Q ₁ (Quantifier)	Q ₂ (Qualifier)
Leucomalachite Green	11.8	330.9	316.2	239.1
Crystal Violet	14.8	372.2	356.3	251.1

Isothiazolinones

MN Appl. No. 128440


Chromatographic conditions

	Column: EC 50/4.6 NUCLEODUR® π ² , 3 μm
	MN REF: 760633.46
	Eluent A: 0.1 % formic acid in water
	Eluent B: acetonitrile – methanol (80:20, v/v)
	Gradient: in 0.5 min from 0 % to 5.0 % B, in 2.0 min from 5 % to 95 % B, hold 95 % B for 0.5 min, in 0.1 min to 0 % B, hold 0 % B for 1.9 min
	Flow rate: 1.3 mL/min
	Temperature: 30 °C
	Detection: MS, MRM
	Injection: 5 μL

Concentration


100 ng/mL in water for each analyte

Chromatograms



MRM transitions


Analyte	RT [min]	[M+H] ⁺	Q ₁ (Quantifier)	Q ₂ (Qualifier)
2-Methyl-4-isothiazolin-3-one (MIT)	1.43	116.1	53.1	45.0
5-Chloro-2-methyl-4-isothiazolin-3-one (CMIT)	1.81	149.9	87.0	135.0
1,2-Benzisothiazol-3(2H)-one (BIT)	1.91	152.0	109.0	134.1
2-Butyl-1,2-benzisothiazolin-3-one (BBIT)	2.52	208.1	152.0	109.0
2-Octyl-4-isothiazolin-3-one (OIT)	2.67	214.2	102.0	43.1
4,5-Dichloro-2-octylisothiazolinone (DCOIT)	3.04	282.1	169.8	73.0



Drug analytes

MN Appl. No. 128380

Chromatographic conditions


Column: EC 50/4.6 NUCLEODUR® π², 3 μm
MN REF: 760633.46
Eluent A: 0.1 % formic acid in water
Eluent B: 0.1 % formic acid in methanol
Gradient: in 4.5 min from 5 % to 90 % B, hold for 0.5 min, in 0.5 min to 5 % B, hold 5 % B for 4.5 min
Flow rate: 1.3 mL/min
Temperature: 30 °C
Detection: MS, MRM
Injection: 5 μL

Concentration



50 ng/mL each in water for each analyte

MRM transitions

Analyte	RT [min]	[M+H] ⁺	Q ₁ (Quantifier)	Q ₂ (Qualifier)
Amphetamine	1.30	136.0	91.1	108.9
3,4-Methylenedioxyamphetamine	1.58	180.0	163.1	105.0
Isoxsuprine	2.58	303.0	285.1	77.1
Fentanyl	3.16	337.0	304.9	105.1
Sibutramine	3.47	280.0	125.0	139.1
Fenbutrazate	3.89	368.2	191.1	91.1

Analyte	RT [min]	[M-H] ⁻	Q ₁ (Quantifier)	Q ₂ (Qualifier)
Hydrochlorothiazide	2.24	295.9	268.7	98.9
Furosemide	4.06	329.0	283.2	255.2

Chromatograms



Nitrosamines

MN Appl. No. 128360

Chromatographic conditions

- Column: EC 100/2 NUCLEODUR® π², 3 μm
- MN REF: 760636.20
- Eluent A: 0.1 % formic acid in water
- Eluent B: 0.1 % formic acid in methanol
- Gradient: hold 2.5 % B for 1.0 min, in 1 min from 2.5 % to 50 % B, hold for 1.0 min, in 0.5 min to 97.5 % B, hold for 2.5 min, in 0.1 min to 2.5 % B, hold for 3.0 min
- Flow rate: 0.3 mL/min
- Temperature: 40 °C
- Detection: MS, MRM
- Injection: 20 μL


Concentration

100 ng/mL in water for each analyte

MRM transitions

Analyte	RT [min]	[M+H] ⁺	Q ₁ (Quantifier)	Q ₂ (Qualifier)
N-Nitroso-diethanolamine	1.64	135.1	74.0	104.1
N-Nitroso-dimethylamine	2.61	75.1	58.0	
N-Nitroso-morpholine	3.22	117.0	45.0	86.0
N-Nitroso-diethylamine	3.94	103.1	75.1	47.0

Chromatogram



Analyte	RT [min]	[M+H] ⁺	Q ₁ (Quantifier)	Q ₂ (Qualifier)
N-Nitroso-piperidine	4.41	115.1	69.1	41.0
N-Nitroso-dipropylamine	4.77	131.1	89.2	43.0
N-Nitroso-dibutylamine	5.02	159.2	57.1	103.2

THC and its metabolites

MN Appl. No. 128640

Chromatographic conditions

- Column: EC 50/4.6 NUCLEODUR® π², 3 μm
- MN REF: 760633.46
- Eluent A: 0.1 % formic acid in water
- Eluent B: 0.1 % formic acid in acetonitrile
- Gradient: in 2.5 min from 0 % to 90 % B, hold for 0.5 min, in 0.1 min to 0 % B, hold 0 % B for 2.9 min
- Flow rate: 1.3 mL/min
- Temperature: 40 °C
- Detection: MS, SMRM
- Injection: 5 μL


Concentration

50 ng/mL in water for each analyte

MRM transitions

Analyte	RT [min]	[M+H] ⁺	Q ₁ (Quantifier)	Q ₂ (Qualifier)
THC-OH	2.64	331.2	99.1	43.1
THC-COOH	2.64	345.2	327.2	299.2
THC	3.06	315.2	193.1	123.1

Chromatograms




Q₁ (Quantifier) — blue line
Q₂ (Qualifier) — red line


Pesticides

MN Appl. No. 128650

Chromatographic conditions

 Column: EC 50/4.6 NUCLEODUR® π², 3 μm
 MN REF: 760633.46
 Eluent A: 0.1 % formic acid in water
 Eluent B: 0.1 % formic acid in methanol
 Gradient: in 5 min from 5 % to 100 % B, hold for 1.0 min. in 0.1 min to 5 % B, hold 5 % B for 3.9 min
 Flow rate: 1.3 mL/min
 Temperature: 30 °C
 Detection: MS, SMRM
 Injection: 5 μL


Fast analysis


 LC/MS analysis of more than 170 pesticides in less than 6 minutes on NUCLEODUR® π², 3 μm.

Concentration

50 ng/mL in water for each analyte

Chromatogram



MRM transitions

Analyte	RT [min]
Cyromazine	0.72
Methamidophos	1.27
Aminocarb	1.40
Acephate	1.55

Analyte	RT [min]
Pymetrozine	1.79
Omethoate	1.88
Mexacarbate	2.20
Thiabendazole	2.26

Continued on page 9

Continued from page 8

Analyte	RT [min]
Fuberidazole	2.27
Aldicarb sulfone	2.35
Nitenpyram	2.35
Monocrotophos	2.55
Flonicamid	2.58
Oxamyl	2.58
Methomyl	2.67
Ethirimol	2.70
Dicrotophos	2.92
Fenuron	3.10
3-Hydroxycarbofuran	3.13
Mevinphos (mix of isomers)	3.20
Thiamethoxam	3.25
Clothianidin	3.26
Prometon	3.26
Dioxacarb	3.30
Secbumeton	3.30
Terbumeton	3.30
Vamidothion	3.34
Dimethoate	3.39
Pirimicarb	3.41
Simetryn	3.43
Imazalil	3.55
Spiroxamine (mix of isomers)	3.62
Aldicarb	3.65
Butocarboxim	3.66
Carbetamide	3.67
Ametryn	3.75
Fenpropimorph	3.75
Imidacloprid	3.79
Propoxur	3.94
Bendiocarb	3.96
Acetamiprid	4.00
Thidiazuron	4.00
Methoprotryne	4.01
Prometryn	4.03
Tebuthiuron	4.04
Carbofuran	4.06
Fluometuron	4.06
Terbutryn	4.09
Carbaryl	4.12
Propham	4.13
Isoprocarb	4.15
Pyracarbolid	4.16
Metribuzin	4.17

Analyte	RT [min]
Monolinuron	4.18
Chlorotoluron	4.19
Isoproturon	4.26
Ethiofencarb	4.27
Cycluron	4.33
Forchlorfenuron	4.34
Hydramethylnon	4.35
Pyrimethanil	4.35
Carboxin	4.36
Fenobucarb	4.39
Oxadixyl	4.39
Isocarbophos	4.42
Flutriafol	4.43
Siduron	4.45
Diuron	4.46
Thiacloprid	4.48
Fludioxinil	4.53
Promecarb	4.53
Spirotetramat	4.54
Halofenozide	4.55
Tricyclazole	4.55
Diethofencarb	4.56
Emamectin B1a benzoate	4.57
Spinosad (Spinosyn A)	4.57
Linuron	4.60
Metobromuron	4.60
Iprovalicarb (mix of isomers)	4.62
Ethiprole	4.63
Metalaxyl	4.63
Paclobutrazol	4.63
Spinetoram	4.64
Spinosad (Spinosyn D)	4.64
Triadimenol	4.64
Butoxycarboxim	4.65
Methabenzthiazuron	4.65
Flutolanil	4.66
Mesotrione	4.66
Methiocarb	4.68
Fipronil	4.70
Cyprodinil	4.74
Cyproconazole (mix of isomers)	4.75
Methoxyfenozide	4.75
Bupirimate	4.76
Fenamidone	4.76
Mepronil	4.77

Continued on page 10

Continued from page 9

Analyte	RT [min]
Chlorantraniliprole	4.78
Ethofumesate	4.82
Fenhexamid	4.82
Furalaxyl	4.82
Triadimefon	4.82
Neburon	4.84
Flubendiamide	4.86
Tetraconazole	4.88
Chloroxuron	4.90
Clethodim (mix of isomers)	4.90
Bifenazate	4.91
Myclobutanil	4.91
Triticonazole	4.93
Mepanipyrim	4.94
Carfentrazone-ethyl	4.99
Diflubenzuron	4.99
Hexaconazole	4.99
Tebuconazole	5.00
Triflumuron	5.04
Diniconazole	5.05
Flufenacet	5.06
Dimoxystrobin	5.10
Fenoxycarb	5.11
Penconazole	5.11
Picoxystrobin	5.12
Flusilazole	5.13
Metconazole	5.13
Butafenacil	5.15
Triflumizole	5.15
Azoxystrobin	5.16
Epoxiconazole	5.16
Mefenacet	5.20
Etaconazole (mix of isomers)	5.21
Kresoxim-methyl	5.24
Diclobutrazol	5.25
Amitraz	5.26
Benalaxyl	5.26

Analyte	RT [min]
Ipconazole (mix of isomers)	5.26
Buprofezin	5.27
Thiobencarb	5.27
Pencycuron	5.29
Bromucanazole (mix of isomers)	5.30
Tebufenpyrad	5.30
Trifloxystrobin	5.31
Cyazofamid	5.32
Fenbuconazole	5.32
Fluoxastrobin	5.32
Prochloraz	5.38
Acibenzolar-S-methyl	5.39
Propiconazole (mix of isomers)	5.40
Clofentezine	5.42
Etoazole	5.44
Flufenoxuron	5.45
Indoxacarb	5.45
Benzoximate	5.46
Tebufenozide	5.47
Piperonyl butoxide	5.49
Propargite	5.51
Furathiocarb	5.52
Spirodiclofen	5.52
Pyraclostrobin	5.53
Pyriproxyfen	5.55
Difenoconazole (mix of isomers)	5.64
Eprinomectin	5.69
Fenpyroximate	5.71
Mandipropamid	5.71
Hexythiazox	5.72
Temephos	5.73
Fenazaquin	5.78
Quinoxifen	5.78
Doramectin	5.84
Pyridaben	5.91
Rotenone	6.00


Chromatogram on page 8



Sulfa drugs

MN Appl. No. 128660

Chromatographic conditions

	Column:	EC 50/4.6 NUCLEODUR® π ² , 3 μm
	MN REF:	760633.46
	Eluent A:	0.1 % formic acid in water
	Eluent B:	0.1 % formic acid in methanol
	Gradient:	in 1.0 min from 5 % to 20 % B, in 4.0 min to 95 % B, hold 95 % B for 0.5 min, in 0.1 min to 5 % B, hold 5 % B for 4.4 min
	Flow rate:	1.3 mL/min
	Temperature:	50 °C
	Detection:	MS, MRM
	Injection:	5 μL

Concentration



100 ng/mL in water for each analyte

MRM transitions

Analyte	RT [min]	[M+H] ⁺	Q ₁ (Quantifier)	Q ₂ (Qualifier)
Sulfafiazole	2.44	256.2	156.2	92.1
Sulfapyridine	2.53	250.2	156.1	92.0
Sulfamethoxypyridazine	2.98	281.2	156.1	92.2
Sulfamonomethoxine	3.15	281.2	156.1	92.2
Sulfachlorpyridazine	3.31	285.2	156.1	92.1

Analyte	RT [min]	[M+H] ⁺	Q ₁ (Quantifier)	Q ₂ (Qualifier)
Sulfacetamide	1.87	215.2	156.2	92.1
Sulfadiazine	2.36	251.2	156.1	92.1
Sulfamerazine	2.75	265.1	156.1	92.1
Sulfadimidine	3.11	279.2	185.9	65.0
Sulfamethoxazole	3.11	254.2	156.1	92.1
Sulfadoxine	3.36	311.1	156.1	92.1
Sulfadimethoxine	3.72	311.1	156.1	92.1
Sulfaquinoxaline	3.84	301.2	156.1	92.1

Chromatograms



NUCLEODUR® π²

Ordering information

Length	50 mm	75 mm	100 mm	125 mm	150 mm	250 mm	EC guard columns*
NUCLEODUR® π², 3 μm							
EC columns (pack of 1)							
2 mm ID	760633.20	760635.20	760636.20	760637.20	760638.20	760639.20	761811.20 (4/2)
3 mm ID	760633.30	760635.30	760636.30	760637.30	760638.30	760639.30	761811.30 (4/3)
4 mm ID	760633.40	760635.40	760636.40	760637.40	760638.40	760639.40	761811.30 (4/3)
4.6 mm ID	760633.46	760635.46	760636.46	760637.46	760638.46	760639.46	761811.30 (4/3)
NUCLEODUR® π², 5 μm							
EC columns (pack of 1)							
2 mm ID	760620.20	760621.20	760622.20	760623.20	760624.20	760625.20	761810.20 (4/2)
3 mm ID	760620.30	760621.30	760622.30	760623.30	760624.30	760625.30	761810.30 (4/3)
4 mm ID	760620.40	760621.40	760622.40	760623.40	760624.40	760625.40	761810.30 (4/3)
4.6 mm ID	760620.46	760621.46	760622.46	760623.46	760624.46	760625.46	761810.30 (4/3)



Selection of guard columns

Guard columns for EC columns with ID	2 mm	3 mm	4 mm	4.6 mm
*EC guard columns (pack of 3)	4 x 2 mm	4 x 3 mm	4 x 3 mm	4 x 3 mm
Required guard column holder (Column protection system)	718966	718966	718966	718966



Guard column holder (Column protection system)



EC guard columns

Registered trademarks

NUCLEODUR® MACHERY-NAGEL GmbH & Co. KG (Germany)
Pinnacle® Restek Corporation (USA)

Distributed By



Greyhound Chromatography and Allied Chemicals
6 Kelvin Park
Birkenhead
Merseyside, CH41 1LT

Tel: 0151 649 4000 Fax: 0151 649 4001
Email: info@greyhoundchrom.com
Web: <https://www.greyhoundchrom.com>

www.mn-net.com

MACHERY-NAGEL

