



Residual Solvent Standard Mixtures for Cannabis Testing

Residual Solvent Mix 1 (16 Components)

1,000 µg/mL in N,N-Dimethylacetamide

Part No.	Price
M-CANRSMIX1AQ4-1ML	\$80.00

Component	CAS No.	Part No.
Acetone	67-64-1	N-11014
Acetonitrile	75-05-8	N-11018
Benzene	71-43-2	N-11149
Carbon tetrachloride	56-23-5	N-11407
Chloroform	67-66-3	N-11447
Ethyl alcohol	64-17-5	N-11885
n-Heptane	142-82-5	N-12543
n-Hexane	110-54-3	N-12550
Isopropyl alcohol	67-63-0	N-12267
Methanol	67-56-1	N-12395
n-Pentane	109-66-0	N-12591
Tetrahydrofuran	109-99-9	N-13545
Toluene	108-88-3	N-13580
m-Xylene	108-38-3	N-12345
o-Xylene	95-47-6	N-12699
p-Xylene	106-42-3	N-12806

Canadian Cannabis Residual Solvent Mixture 1 (26 Components)

5,000 µg/mL in N,N-Dimethylacetamide

Part No.	Price
M-CANRSMIX1AQ7-1ML	\$175.00

Component	CAS No.	Part No.
Acetic acid	64-19-7	N-11011
Acetone	67-64-1	N-11014
n-Amyl alcohol	71-41-0	N-12509
Anisole	100-66-3	N-11080
2-Butanone	78-93-3	N-10297
n-Butyl acetate	123-86-4	N-12512
sec-Butyl alcohol	78-92-2	N-13204
n-Butyl alcohol	71-36-3	N-12514
tert-Butyl methyl ether	1634-04-4	N-13528
Ethyl acetate	141-78-6	N-11881
Ethyl alcohol	64-17-5	N-11885
Ethyl ether	60-29-7	N-11897
Ethyl formate	109-94-4	N-11898
Formic acid	64-18-6	N-12016
n-Heptane	142-82-5	N-12543
Isoamyl alcohol	123-51-3	N-12225
Isobutyl acetate	110-19-0	N-12233
Isobutyl alcohol	78-83-1	N-12234
Isopropyl acetate	108-21-4	N-12223
Isopropyl alcohol	67-63-0	N-12267
Methyl acetate	79-20-9	N-12411
Methyl sulfoxide	67-68-5	N-12457
n-Pentane	109-66-0	N-12591
Propyl acetate	109-60-4	N-13120
Propyl alcohol	71-23-8	N-13121
Triethylamine	121-44-8	N-13674

California Residual Solvent Mixture 1 (20 Components)

Varied Concentration in N,N-Dimethylacetamide

Part No.	Price
M-CARSMIX1AQ99-1ML	\$200.00

Component	CAS No.	Concentration	Part No.
Acetone	67-64-1	28,000 µg/mL	N-11014
Acetonitrile	75-05-8	4,920 µg/mL	N-11018
Benzene	71-43-2	24 µg/mL	N-11149
Chloroform	67-66-3	720 µg/mL	N-11447
1,2-Dichloroethane	107-06-2	60 µg/mL	N-10154
N,N-Dimethylacetamide	127-19-5	937,000 µg/mL	N-12625
Ethyl acetate	141-78-6	28,000 µg/mL	N-11881
Ethyl alcohol	64-17-5	28,000 µg/mL	N-11885
Ethyl ether	60-29-7	28,000 µg/mL	N-11897
n-Heptane	142-82-5	28,000 µg/mL	N-12543
n-Hexane	110-54-3	3,480 µg/mL	N-12550
Isopropyl alcohol	67-63-0	28,000 µg/mL	N-12267
Methanol	67-56-1	36,000 µg/mL	N-12395
Methylene chloride	75-09-2	7,200 µg/mL	N-12471
n-Pentane	109-66-0	28,000 µg/mL	N-12591
Toluene	108-88-3	10,680 µg/mL	N-13580
Trichloroethene	79-01-6	960 µg/mL	N-13654
m-Xylene	108-38-3	26,040 µg/mL	N-12345
o-Xylene	95-47-6	26,040 µg/mL	N-12699
p-Xylene	106-42-3	26,040 µg/mL	N-12806

Oregon Residual Solvent Mixture 1 (20 Components)

Varied Concentration in N,N-Dimethylacetamide

Part No.	Price
M-CARSMIX1AQ99-1ML	\$200.00

Component	CAS No.	Concentration	Part No.
Acetone	67-64-1	28000 ug/mL	N-11014
Acetonitrile	75-05-8	4920 ug/mL	N-11018
Benzene	71-43-2	24 ug/mL	N-11149
Chloroform	67-66-3	720 ug/mL	N-11447
1,2-Dichloroethane	107-06-2	60 ug/mL	N-10154
N,N-Dimethylacetamide	127-19-5	937000 ug/mL	N-12625
Ethyl acetate	141-78-6	28000 ug/mL	N-11881
Ethyl alcohol	64-17-5	28000 ug/mL	N-11885
Ethyl ether	60-29-7	28000 ug/mL	N-11897
n-Heptane	142-82-5	28000 ug/mL	N-12543
n-Hexane	110-54-3	3480 ug/mL	N-12550
Isopropyl alcohol	67-63-0	28000 ug/mL	N-12267
Methanol	67-56-1	36000 ug/mL	N-12395
Methylene chloride	75-09-2	7200 ug/mL	N-12471
n-Pentane	109-66-0	28000 ug/mL	N-12591
Toluene	108-88-3	10680 ug/mL	N-13580
Trichloroethene	79-01-6	960 ug/mL	N-13654
m-Xylene	108-38-3	26040 ug/mL	N-12345
o-Xylene	95-47-6	26040 ug/mL	N-12699
p-Xylene	106-42-3	26040 ug/mL	N-12806

Distributed By:



Greyhound Chromatography and Allied Chemicals
6 Kelvin Park
Birkenhead
Merseyside, CH41 1LT

Tel: 0151 649 4000 Fax: 0151 649 4001
Email: info@greyhoundchrom.com
Web: <https://www.greyhoundchrom.com>