

WELLINGTON Reporter

PRODUCT UPDATES FROM WELLINGTON LABORATORIES

February 20, 2013

NEW PRODUCTS

Native and Mass-labelled Hexafluoropropylene Oxide Dimer Acid

Hexafluoropropylene oxide (HFPO) is a well-known versatile synthetic building block in the manufacturing of fluoropolymers (such as perfluoroalkoxy plastics) as well as a number of poly- and per-fluorinated intermediates. Although it is used to produce a vast number of commercial products, it's reactivity makes its survival in the environment unlikely. However, HFPO can react to form a stable dimer acid during oligomerization, or other manufacturing processes, which could lead to its detection in environmental samples. The presence of this HFPO dimer acid (HFPO-DA) in the environment could be due to residual leaching from commercial products or direct release during the manufacturing processes.

For this reason, Wellington has synthesized a native and mass-labelled ($^{13}C_3$) hexafluoropropylene oxide dimer acid reference standard, HFPO-DA and M3HFPO-DA respectively, to aid researchers in their quantification of this potential environmental contaminant.

HFPO-DA

M3HFPO-DA

Catalogue Number	Product (methanol)	Qty/Conc
HFPO-DA	2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)propanoic acid	1.2 ml 50 µg/ml
МЗНГРО-DA	2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)-13C3-propanoic acid	1.2 ml 50 µg/ml





Greyhound Chromatography & Allied Chemicals 6 Kelvin Park, Birkenhead, Merseyside, CH41 1LT, U.K. Tel: (+44)-0-151-649-4000 Fax: (+44)-0-151-649-4001 E-mail: info@ greyhoundchrom.com Web: www.greyhoundchrom.com

Certified System

_WELLINGTON LABORATORIES