




NEW PRODUCT

Sodium perfluoro-1-propanesulfonate **L-PFPrS**

Aqueous Film Forming Foams (AFFF) are a group of fire-fighting foams containing perfluoroalkyl substances (PFAS) that have been utilized for fire extinguishment and vapour suppression by facilities handling large quantities of flammable hydrocarbon liquids (e.g. airports, storage terminals, refineries, etc.) since the 1970s. It has been documented that perfluoro-1-propanesulfonate can be present in some AFFF formulations at levels as high as 3.5% of the total perfluoroalkanesulfonic acid (PFSA) concentration. In fact, recent detection of ultra-short PFASs in ground water has elevated scientific interest in current environmental levels, steady-state human exposure, and the environmental fate of these compounds. For this reason, **Wellington** has prepared an analytical reference standard for sodium perfluoro-1-propanesulfonate (**L-PFPrS**) to aid researchers in both method development and quantification.



Sodium perfluoro-1-propanesulfonate

	Catalogue Number	Product (methanol)	Qty/Conc	
NEW	L-PFPrS	Sodium perfluoro-1-propanesulfonate	1.2 ml	50 µg/ml
	L-PFBS	Potassium perfluoro-1-butanesulfonate	1.2 ml	50 µg/ml
	L-PFPeS	Sodium perfluoro-1-pentanesulfonate	1.2 ml	50 µg/ml
	L-PFHxS	Sodium perfluoro-1-hexanesulfonate	1.2 ml	50 µg/ml
	L-PFHpS	Sodium perfluoro-1-heptanesulfonate	1.2 ml	50 µg/ml
	L-PFOS	Sodium perfluoro-1-octanesulfonate	1.2 ml	50 µg/ml
	L-PFOSK	Potassium perfluoro-1-octanesulfonate	1.2 ml	50 µg/ml
	L-PFNS	Sodium perfluoro-1-nonanesulfonate	1.2 ml	50 µg/ml
	L-PFDS	Sodium perfluoro-1-decanesulfonate	1.2 ml	50 µg/ml
	L-PFDoS	Sodium perfluoro-1-dodecanesulfonate	1.2 ml	50 µg/ml



ISO 9001

Distributed Throughout Europe and Middle East By-

Greyhound Chromatography & Allied Chemicals
6 Kelvin Park, Birkenhead,
Merseyside, CH41 1LT, U.K.
Tel: (+44)-0-151-649-4000 Fax: (+44)-0-151-649-4001
E-mail: info@greyhoundchrom.com
Web: www.greyhoundchrom.com

