ADAMDU Software



Available as a download

Streamlines data collection from multiple laboratory balances

Features

- Monitor and collect data readings for up to eight balances or scales simultaneously
- · Include date and time with readings
- · Graph readings in real time
- Export readings in several file formats (text, CSV, HTML, XML)
- Export directly to Microsoft® Excel or Microsoft® Word

Microsoft Word, Microsoft Excel, and Microsoft are all registered trademark of Microsoft Corporation in the United States and other countries

Part Number	Class
600002028	Adam DU

Anti-Vibration Table



Sturdy base helps minimise movement and reduce fluctuations in measurement

Features

- · Scratch-resistant granite slab
- Two-table construction keeps the granite block in the middle and reduces interference
- Large working surface of 400x450mm
- Area surrounding the slab for placing samples
- · Lightweight tubular aluminum frame
- Adjustable levelling feet for uneven surfaces
- Polished finish

Part Number	Class
104008036	Anti-Vibration Table