

DDK Density Determination Kits

Calculate densities for solid and liquid samples



Features

- Glass beaker
- Suspended weighing cradle to replace the standard pan
- Independent support for beaker
- Thermometer with clip
- Glass sinker of known volume for the determination of liquid density

Part Number	Class
1060012713	DDK Density Determination Kit (90mm pan size)
1060012714	DDK Density Determination Kit (120 and 160mm pan size)
1060014262	Density kit for Equinox and Solis

ATP Printer

Print results quickly and easily

Features

- Thermal printer
- Dimensions: 114 mm x 184 mm x 90 mm
- Print size 24 x 24 dot
- Column size 384 dots per line
- 57.5mm paper width



Part Number	Class
1120011156	ATP Printer
3126011281	ATP Printer Paper (10 Rolls)
3126001263	ATP Printer Paper